


Programming Guide

SK305 *Aigoual* Linear TEC Driver

SK-Series Modules



General Information

Important Notice

Information in this document is subject to change without notice.
Copyright © SISYPH, 2025. All rights reserved.

Signals and Systems for Physics
BP90406
16 place saint-Georges
F31000 Toulouse France
Phone (+33) 781 547 391
www.sisyph.com

Scope

This document describes operating the SK305 Linear TEC Driver module over the serial interface.

PRELIMINARY INFORMATION

Contents

General Information	2
Important Notice	2
Scope	2
1 Introduction	4
1.1 Power-on Configuration	4
1.2 Buffers	4
1.3 Command syntax	4
1.4 Examples	5
2 List of Commands	6
2.1 Instrument Settings commands	8
2.2 Instrument Configuration commands	14
2.3 Instrument Monitoring commands	21
2.4 Status Reporting commands	26
2.5 Interface commands	39
2.6 Memory commands	48
3 Status Model	50
3.1 Master Summary Status (MSTS)	51
3.2 Master Summary Enable (MSTE)	51
3.3 Event Status (EVTS)	52
3.4 Event Enable (EVTE)	52
3.5 Instrument Status (INSS)	53
3.6 Instrument Enable (INSE)	53
3.7 Instrument Condition (INSC)	53
3.8 Overload Status (OVLS)	54
3.9 Overload Enable (OVLE)	54
3.10 Overload Condition (OVLG)	54
3.11 Communication Status (COMS)	55
3.12 Communication Enable (COME)	55
3.13 Last Command Error (LCMD)	56
3.14 Last Execution Error (LEXE)	56
3.15 Last Instrument Error (LINS)	56
3.16 Last User Request (LURQ)	56
4 Index of commands	57
5 Document Revision History	58
5.1 Version Number	58
5.2 Revision History	58

1 Introduction

Remote operation of the SK305 is through a simple command language documented in this chapter. Both set and query forms of most commands are supported, allowing the user complete control of the module from a remote computer.

1.1 Power-on Configuration

The settings for serial interface are 9600 baud with no parity and no hardware flow control, and local echo disabled (CONS 0).

Most of the instrument settings are stored in non-volatile memory and can be retrieved using the appropriate commands. At power-on the instrument returns to the state noted in the command descriptions. Reset values (*RST command) of parameters are shown in **boldface**.

1.2 Buffers

The instrument stores incoming bytes from the host interface in a 128-byte input buffer. Characters accumulate in the input buffer until a command terminator (either <CR> or <LF>) is received, at which point the message is parsed and executed. Query responses from the instrument are sent when they are ready without any flow control nor output buffering. The input buffer is automatically flushed upon detecting an overflow, and an error is recorded in the EVTS status register.

1.3 Command syntax

The four letter mnemonic (shown in CAPS) in each command sequence specifies the command. The rest of the sequence consists of parameters. The command parser accepts only uppercase mnemonics.

Commands may take either set or query form, depending on whether the ? character follows the mnemonic. *Set only* commands are listed without the ?, *query only* commands show the ? after the mnemonic, and *optionally query* commands are marked with a (?). Parameters shown in { } and [] are not always required. Parameters in { } are only required to set a value, and should be omitted for queries. Parameters in [] are optional in both set and query commands. Parameters listed without any surrounding characters are always required. Do not send () or { } or [] as part of the command. Multiple parameters are separated by commas. Multiple commands may be sent on one command line by separating them with semicolons ; so long as the input buffer does not overflow. Commands are terminated by either <CR> or <LF> characters. Null commands and whitespace are ignored. Execution of the command does not begin until the command terminator is received.

The following table summarizes the notation used in the command descriptions:

Symbol	Definition
<i>b</i>	boolean
<i>i, m, n</i>	unsigned integers
<i>u, v</i>	signed integers
(?)	required for queries; illegal for set commands.
<i>p</i>	parameter always required.
{ <i>p</i> }	required parameter for set commands; illegal for queries.
[<i>p</i>]	optional parameter for both set and query forms.

1.4 Examples

Each command is provided with a simple example illustrating its usage. In these examples, all data sent by the host computer to the instrument are set as *straight teletype font*, while responses received the host computer from the instrument are set as *slanted teletype font*. The usage examples vary with respect to set/query, optional parameters, and token formats. These examples are not exhaustive, but are intended to provide a convenient starting point for user programming.

PRELIMINARY INFORMATION

2 List of Commands

This section provides syntax and operational descriptions for remote commands.

2.1 Instrument Settings commands	8
MANS (Manual Current Control)	8
ILMP (Positive Current Limit)	9
ILMN (Negative Current Limit)	10
VTHP (Positive Voltage Threshold)	11
VTHN (Negative Voltage Threshold)	12
FFWG (Feedforward Gain)	13
2.2 Instrument Configuration commands	14
MANE (Manual Control State)	14
EXTE (External Control State)	15
FFWE (External Current Control State)	16
TECE (TEC Output State)	17
ITPO (Trip-off upon Current Limiting)	18
VTPO (Trip-off upon Overvoltage Detection)	19
MONS (MONO Output Selection)	20
2.3 Instrument Monitoring commands	21
RMON (Remote Monitoring)	21
TDIE (Die Temperature)	22
STMS (Streamed Channels Selection)	23
STME (Data Streaming Enable)	24
STMN (Number of Streamed Measurements)	25
2.4 Status Reporting commands	26
*CLS (Clear Status Registers)	26
MSTS (Master Summary Status)	27
MSTE (Master Summary Enable)	28
EVTS (Event Status)	29
EVTE (Event Enable)	30
COMS (Communications Status)	31
COME (Communications Enable)	32
OVLS (Overload Status)	33
OVLE (Overload Enable)	34
OVLC (Overload Condition)	35
INSS (Instrument Status)	36
INSE (Instrument Enable)	37
INSC (Instrument Condition)	38
2.5 Interface commands	39
*RST (Reset)	39
*OPC (Operation Complete)	40
CONS (Console Mode)	41
*IDN (Identify)	42
LINS (Last Instrument Error Status)	43

LURQ (Last User Request Status)	44
LCMD (Last Command Error Status)	45
LEXE (Last Execution Error Status)	46
TERM (Response Termination)	47
2.6 Memory commands	48
*RCL (Recall Settings)	48
*SAV (Save Current Settings)	49

PRELIMINARY INFORMATION

2.1 Instrument Settings commands

The Instrument Settings commands provide control of the instrument's physical parameters.

MANS (Manual Current Control)

Group	Instrument Settings commands
Action	Set/Query
Syntax	MANS(?) <i>{u}</i>
Description	Set (query) the setpoint current <i>{to u}</i> , in mA, with a resolution of 12-bit.
Allowed range	$-1000 \leq u \leq 1000$
Power-on value	Restored from non-volatile memory (<i>cf</i> *SAV).
Reset value	0
Example	MANS 500; MANS? <i>500</i>
Related commands	MANE.

PRELIMINARY INFORMATION

ILMP (Positive Current Limit)

Group	Instrument Settings commands
Action	Set/Query
Syntax	ILMP(?) <i>{u}</i>
Description	<p>Set (query) the positive current limit <i>{to u}</i>, in mA, with a resolution of 12-bit.</p> <p>When the output current demand is above the positive limit, the internal limiting circuitry is engaged to protect the load. In this operation, the output current never exceeds the positive limit, whatever the combination of the control inputs. Such an event can also switch the current source off, this latter functionality being controlled through the ITPO command.</p>
Allowed range	$0 \leq u \leq 1000$
Power-on value	Restored from non-volatile memory (cf *SAV).
Reset value	1000
Example	<pre>ILMP 500; ILMP? 500</pre>
Related commands	ITPO, ILMN.

ILMN (Negative Current Limit)

Group	Instrument Settings commands
Action	Set/Query
Syntax	ILMN(?) <i>{u}</i>
Description	<p>Set (query) the negative current limit <i>{to u}</i>, in mA, with a resolution of 12-bit.</p> <p>When the output current demand is below the negative limit, the internal limiting circuitry is engaged to protect the load. In this operation, the output current never exceeds the negative limit, whatever the combination of the control inputs. Such an event can also switch the current source off, this latter functionality being controlled through the ITPO command.</p>
Allowed range	$-1000 \leq u \leq 0$
Power-on value	Restored from non-volatile memory (of *SAV).
Reset value	-1000
Example	<pre>ILMN -500; ILMN?</pre> <p><i>-500</i></p>
Related commands	ITPO, ILMP.

VTHP (Positive Voltage Threshold)

Group	Instrument Settings commands
Action	Set/Query
Syntax	VTHP(?) { <i>u</i> }
Description	<p>Set (query) the positive voltage threshold {to <i>u</i>}, in mV, with a resolution of 12-bit.</p> <p>The compliance voltage of the current source is intrinsically limited about ± 4.5 V. Nevertheless, the ILMP command can be used to set the positive limit to a lower value. When the output voltage is above the positive threshold, the corresponding flag is raised in the OVLS status register. This event can also switch the current source off, this latter functionality being controlled through the VTPO command.</p>
Allowed range	$0 \leq u \leq 5000$
Power-on value	Restored from non-volatile memory (cf *SAV).
Reset value	5000
Example	VTHP 2500; VTHP? 2500
Related commands	VTPO, VTHN.

VTHN (Negative Voltage Threshold)

Group	Instrument Settings commands
Action	Set/Query
Syntax	VTHN(?) { <i>u</i> }
Description	<p>Set (query) the negative voltage threshold {to <i>u</i>}, in mV, with a resolution of 12-bit.</p> <p>The compliance voltage of the current source is intrinsically limited about ± 4.5 V. Nevertheless, the ILMN command can be used to set the negative limit to an upper value. When the output voltage is below the negative threshold, the corresponding flag is raised in the OVLS status register. This event can also switch the current source off, this latter functionality being controlled through the VTPO command.</p>
Allowed range	$-5000 \leq u \leq 0$
Power-on value	Restored from non-volatile memory (cf *SAV).
Reset value	-5000
Example	VTHN -2500; VTHN? <i>-2500</i>
Related commands	VTPO, VTHP.

PRELIMINARY INFORMATION

FFWG (Feedforward Gain)

Group	Instrument Settings commands
Action	Set/Query
Syntax	FFWG(?) { <i>u</i> }
Description	<p>Set (query) the transmission coefficient of the feedforward input {to <i>u</i>}, in ‰, with a resolution of 12-bit.</p> <p>This command is used to adjust both weight and polarity of the feed-forward input voltage. Indeed, the command FFWG 1000 will set the adjustable transmission gain to +1 V/V, while FFWG 0 will provide maximal attenuation. In order to reverse the transmission polarity, the command is invoked with negative argument, for instance, FFWG -1000 will just multiply the FFWI signal by the scalar -1.</p>
Allowed range	$-1000 \leq u \leq 1000$
Power-on value	Restored from non-volatile memory (of *SAV).
Reset value	0
Example	FFWG 100; FFWG? <i>100</i>
Related commands	FFWE.

PRELIMINARY INFORMATION

2.2 Instrument Configuration commands

The Instrument Configuration commands provide control of the instrument's physical functionalities.

MANE (Manual Control State)

Group	Configuration commands
Action	Set/Query
Syntax	MANE(?) { <i>b</i> }
Description	Set (query) the manual control's state {to <i>b</i> }. When <i>b</i> = 1 (resp. 0), the manual current control functionality is enabled (resp. enabled). The MANS command should be invoked prior enabling the current source output.
Allowed range	$b \in \{0, 1\}$
Power-on value	Restored from non-volatile memory (cf. *SAV).
Reset value (*RST)	1
Example	MANE 1; MANE? <i>1</i>
Related commands	MANS

EXTE (External Control State)

Group	Configuration commands
Action	Set/Query
Syntax	EXTE(?) { <i>b</i> }
Description	Set (query) the external control's state {to <i>b</i> }. When <i>b</i> = 1 (resp. 0), the external current control input is enabled (resp. enabled).
Allowed range	$b \in \{0, 1\}$
Power-on value	Restored from non-volatile memory (cf *SAV).
Reset value (*RST)	0
Example	EXTE 1; EXTE? <i>1</i>
Related commands	

PRELIMINARY INFORMATION

FFWE (Feedforward Control State)

Group	Configuration commands
Action	Set/Query
Syntax	FFWE(?) { <i>b</i> }
Description	Set (query) the feedforward control's state {to <i>b</i> }. When <i>b</i> = 1 (resp. 0), the feedforward current control input is enabled (resp. enabled).
Allowed range	$b \in \{0, 1\}$
Power-on value	Restored from non-volatile memory (cf *SAV).
Reset value (*RST)	0
Example	FFWE 1; FFWE? <i>1</i>
Related commands	FFWG

PRELIMINARY INFORMATION

TECE (TEC Output State)

Group	Configuration commands
Action	Set/Query
Syntax	TECE(?) { <i>b</i> }
Description	<p>Set (query) the current source output's state {to <i>b</i>}.</p> <p>When $b = 1$ (resp. 0), the TEC output is enabled (resp. disabled) and the programmed current is delivered from the source to the load. Otherwise, both load and current source terminals are shorted and tied to the ground voltage through relay contacts.</p>
Allowed range	$b \in \{0, 1\}$
Power-on value	Restored from non-volatile memory (cf *SAV).
Reset value (*RST)	0
Example	TECE 1; TECE? <i>1</i>
Related commands	

PRELIMINARY INFORMATION

ITPO (Trip-off upon Current Limiting)

Group	Instrument Configuration commands
Action	Set/Query
Syntax	ITPO(?) { <i>n</i> }
Description	<p>Set (query) the current source to trip-off upon current limiting {to <i>n</i>}.</p> <p>The programmable limiter restricts the range of the control voltage driving the output transconductance amplifier. This limiter prevents excessive amplifier's input voltage using analog circuitry and DAC voltages programmed <i>via</i> the ILMP and ILMN commands. Two events, reporting that the positive (or negative) voltage limiter is engaged, can be used to switch the current source off. For certain critical situations, this functionality can be useful to detect abnormal operating conditions and automatically shut the current source down. The ITPO command is provided to select what event is monitored to trip-off upon current limiting.</p>
Allowed range	<p>$n \in \{0, \dots, 3\}$ where :</p> <ul style="list-style-type: none"> 0 \longleftrightarrow normal operation (current limiting); 1 \longleftrightarrow tripped-off upon positive current limiting; 2 \longleftrightarrow tripped-off upon negative current limiting; 3 \longleftrightarrow tripped-off upon positive and negative current limiting.
Power-on value	Restored from non-volatile memory (<i>cf</i> *SAV).
Reset value	0
Example	ITPO 1
Related commands	ILMP, ILMN.

VTPO (Trip-off upon Overvoltage Detection)

Group	Instrument Configuration commands
Action	Set/Query
Syntax	VTPO(?) { <i>n</i> }
Description	<p>Set (query) the current source to trip-off upon over-voltage detection {to <i>n</i>}.</p> <p>The current source output is continuously monitored to detect excessive compliance voltage. For this purpose, the VTHP and VTHN commands provide programmable threshold levels for the positive and negative detectors. While no specific action (<i>e.g.</i> voltage limiting) is carried on when an out-of-range compliance voltage is detected, two events, reporting that the positive (or negative) output voltage limit is reached, can be used to switch the current source off. For certain critical situations, this functionality can be useful to detect abnormal operating conditions and automatically shut the current source down. The ITPO command is provided to select what event is monitored to trip-off upon over-voltage detection.</p>
Allowed range	<p>$n \in \{0, \dots, \mathbf{3}\}$ where :</p> <ul style="list-style-type: none"> 0 \longleftrightarrow normal operation; 1 \longleftrightarrow tripped-off upon positive detection; 2 \longleftrightarrow tripped-off upon negative detection; 3 \longleftrightarrow tripped-off upon positive and negative detection.
Power-on value	Restored from non-volatile memory (<i>cf</i> *SAV).
Reset value	3
Example	VTPO 1
Related commands	VTHP, VTHN.

MONS (MONO Output Selection)

Group	Configuration commands
Action	Set/Query
Syntax	MONS(?) { <i>n</i> }
Description	Set (query) the source of the monitoring output signal {to <i>n</i> }. This command is used to select which signal is routed to the MONO coaxial connector.
Allowed range	$n \in \{0, \dots, 3\}$ where : 0 ↔ ground voltage (0V); 1 ↔ output current monitor (IMON); 2 ↔ output voltage monitor (VMON); 3 ↔ /STATUS signal.
Power-on value	Restored from non-volatile memory (cf *SAV).
Reset value	0
Example	MONS 1; MONS? 1
Related commands	

PRELIMINARY INFORMATION

2.3 Instrument Monitoring commands

The Instrument Monitoring commands provide the host computer with the last measurements of the instrument's physical parameters.

RMON (Remote Monitoring)

Group	Monitoring commands
Action	Query only
Syntax	RMON? <i>n</i>
Description	Return to the host computer the last measurement of the parameter specified by <i>n</i> . These parameters are periodically sampled (100 ms) by an internal task.
Allowed range	$n \in \{1, 2\}$, where the index <i>n</i> : 1 \longleftrightarrow output current (IMON), in mA; 2 \longleftrightarrow output voltage (VMON), in mV.
Example	RMON? 1 -23
Related commands	

TDIE (Die Temperature)

Group	Monitoring commands
Action	Query only
Syntax	TDIE?
Description	Return the die temperature. TDIE? returns the last measurement of the temperature (in K) of the die provided by the MCU on-chip sensor. The precision is about ± 1 K. This reading can be used to get an approached value of the main printed circuit board's temperature where the MCU is mounted. This measurement is automatically updated every 100 ms.
Example	TDIE? 298
Related commands	

PRELIMINARY INFORMATION

STMS (Streamed Channels Selection)

Group	Monitoring commands
Action	Set/Query
Syntax	STMS(?) { <i>m</i> }
Description	<p>Set (query) the channels selection register {to bit-mask <i>m</i>}.</p> <p>In order to stream the channel <i>i</i> to output, the command must be invoked with $m = 2^i$ as argument. For instance, STMS 6 will be executed to stream the positive and negative peak-voltages of the PI2D command to output. Data are output on a single line where a comma delimiter is used to separate channel data. The rightmost position in the row is occupied by the channel with the lowest weight. Data streaming operation is enabled using the STME command while the number of measurements that will be streamed out is set by the STMN command.</p>
Allowed range	<p>$m \in \{1, \dots, 3\}$, where <i>m</i> can be any combination of :</p> <p>$m_0 = 1 = 2^0 \longleftrightarrow$ output current (IMON), in mA;</p> <p>$m_1 = 2 = 2^1 \longleftrightarrow$ output voltage (VMON), in mV.</p>
Power-on value	Restored from non-volatile memory (cf *SAV).
Reset value	1
Example	<pre>STMN 4; STMS 3 ; STME 1 602,2413 601,2412 602,2413 602,2413</pre>
Related commands	STME, STMN, TDIE.

STME (Data Streaming Enable)

Group	Monitoring commands
Action	Set/Query
Syntax	STME(?) { <i>b</i> }
Description	Set (query) the state of the data streaming functionality {to <i>b</i> }. If <i>b</i> = 0 (resp. 1), data streaming is disabled (resp. enabled). When data streaming is enabled, the number of measurements that will be streamed out is specified using the STMN command.
Allowed range	$b \in \{0, 1\}$.
Power-on value	0
Reset value (*RST)	0
Example	STME 1; STME? <i>1</i>
Related commands	STMS, STMN.

PRELIMINARY INFORMATION

STMN (Number of Streamed Measurements)

Group	Monitoring commands
Action	Set/Query
Syntax	STMN(?) { <i>n</i> }
Description	<p>Set (query) the number of measurements to be streamed out {to <i>n</i>}.</p> <p>When $n = 0$, measurements will be output indefinitely, until the STME command is invoked. If $n > 0$, data streaming will stop once n measurements will have been output. Measurements are streamed out at a constant rate of approximately one measurement per second.</p>
Allowed range	$0 \leq n \leq 10000$.
Reset value	0
Example	<p>STMN 1000; STMN?</p> <p>1000</p>
Related commands	STME, STMS.

PRELIMINARY INFORMATION

2.4 Status Reporting commands

The Status commands query and configure registers associated with status reporting of the instrument.

*CLS (Clear Status Registers)

Group	Status reporting commands
Action	Query only
Syntax	*CLS
Description	Clear immediately all status registers, which are : LEXE, LCMD, LINS, LURQ, INSS, OVLS, COMS and EVTS.
Example	*CLS
Related commands	

PRELIMINARY INFORMATION

MSTS (Master Summary Status)

Group	Status reporting commands
Action	Query only
Syntax	MSTS? [<i>n</i>]
Description	<p>Return the Master Summary Status register [bit-mask <i>n</i>].</p> <p>The execution of the MSTS? query – without the optional bit-mask <i>n</i> – always causes the /STATUS signal to be de-asserted. Note that MSTS? <i>n</i> will not clear /STATUS, even if bit $i \mid n = 2^i$ is the only bit presently causing the /STATUS signal.</p>
Power-on value	0
Example	<pre>MSTS?; MSTS? 128;</pre> <pre>129</pre> <pre>128</pre>
Related commands	MSTE

PRELIMINARY INFORMATION

MSTE (Master Summary Enable)

Group	Status reporting commands
Action	Set/Query
Syntax	MSTE(?) [<i>n</i>] { <i>m</i> }
Description	Set (query) the Master Summary Enable register [bit-mask <i>n</i>] {to bit-mask <i>m</i> }. The set-form command will clear the bits outside the bit-mask.
Power-on value	0
Example	MSTE 128; MSTE? <i>128</i>
Related commands	MSTS

PRELIMINARY INFORMATION

EVTS (Event Status)

Group	Status reporting commands
Action	Query only
Syntax	EVTS? [<i>n</i>]
Description	Read the Event Summary Status register [bit-mask <i>n</i>].
Power-on value	1
Example	EVTS? 4
Related commands	EVTE

PRELIMINARY INFORMATION

EVTE (Event Enable)

Group	Status reporting commands
Action	Set/Query
Syntax	EVTE(?) [<i>n</i>] { <i>m</i> }
Description	Set (query) the Event Summary Enable register [bit-mask <i>n</i>] {to bit-mask <i>m</i> }. The set-form command will clear the bits outside the bit-mask.
Power-on value	0
Example	EVTE 4; EVTE? 4
Related commands	EVTS

PRELIMINARY INFORMATION

COMS (Communications Status)

Group	Status reporting commands
Action	Query only
Syntax	COMS? [<i>n</i>]
Description	Read the Communications Status register [bit-mask <i>n</i>].
Power-on value	0
Example	COMS? <i>0</i>
Related commands	COME

PRELIMINARY INFORMATION

COME (Communications Enable)

Group	Status reporting commands
Action	Set/Query
Syntax	COME(?) [<i>n</i>] { <i>m</i> }
Description	Set (query) the Communications Enable register [bit-mask <i>n</i>] {to bit-mask <i>m</i> }. The set-form command will clear the bits outside the bit-mask.
Power-on value	0
Example	COME 1
Related commands	COMS

PRELIMINARY INFORMATION

OVLS (Overload Status)

Group	Status reporting commands
Action	Query only
Syntax	OVLS? [<i>n</i>]
Description	Read the Overload Status register [bit-mask <i>n</i>].
Power-on value	0
Example	OVLS? 2
Related commands	OVLE, OVLC.

PRELIMINARY INFORMATION

OVLE (Overload Enable)

Group	Status reporting commands
Action	Set/Query
Syntax	OVLE(?) [<i>n</i>] { <i>m</i> }
Description	Set (query) the Overload Enable register [bit-mask <i>n</i>] {to bit-mask <i>m</i> }. The set-form command will clear the bits outside the bit-mask.
Power-on value	0
Example	OVLE 2
Related commands	OVLS, OVLC.

PRELIMINARY INFORMATION

OVLC (Overload Condition)

Group	Status reporting commands
Action	Query only
Syntax	OVLC? [<i>n</i>]
Description	<p>Read the Overload Condition register [bit-mask <i>n</i>].</p> <p>The values of the bits in the OVLC condition register are determined by the current (real-time) condition of the events defined in the OVLS status register.</p> <p>Reading the condition register does not affect the register.</p>
Power-on value	0
Example	<p>OVLC?</p> <p>2</p>
Related commands	OVLS, OVLE.

PRELIMINARY INFORMATION

INSS (Instrument Status)

Group	Status reporting commands
Action	Query only
Syntax	INSS? [<i>n</i>]
Description	Read the Instrument Status register [bit-mask <i>n</i>].
Power-on value	0
Example	INSS? <i>1</i>
Related commands	LINS, INSE, INSC.

PRELIMINARY INFORMATION

INSE (Instrument Enable)

Group	Status reporting commands
Action	Set/Query
Syntax	INSE(?) [<i>n</i>] { <i>m</i> }
Description	Set (query) the Instrument Enable register [bit-mask <i>n</i>] {to bit-mask <i>m</i> }. The set-form command will clear the bits outside the bit-mask.
Power-on value	0
Example	INSE 2
Related commands	LINS, INSS, INSC.

PRELIMINARY INFORMATION

INSC (Instrument Condition)

Group	Status reporting commands
Action	Query only
Syntax	INSC? [<i>n</i>]
Description	<p>Read the Instrument Condition register [bit-mask <i>n</i>].</p> <p>The values of the bits in the INSC condition register are determined by the current (real-time) condition of the events defined in the INSS status register.</p> <p>Reading the condition register does not affect the register.</p>
Power-on value	0
Example	<p>INSC?</p> <p>2</p>
Related commands	LINS, INSE, INSS.

PRELIMINARY INFORMATION

2.5 Interface commands

The Interface commands provide control over the interface between the instrument and the host computer.

*RST (Reset)

Group	Interface commands
Action	Set only
Syntax	*RST
Description	<p>Reset the instrument to its default configuration.</p> <p>When a parameter is affected by the *RST command, its value is reset according to the information given by the Reset value field within the related command section.</p> <p>Whereas status registers are unaffected by *RST, the content of some conditions registers may have been modified upon resetting the instrument.</p>
Example	*RST
Related commands	*RCL, *SAV.

PRELIMINARY INFORMATION

***OPC (Operation Complete)**

Group	Interface commands
Action	Set/Query
Syntax	*OPC(?)
Description	Set the OPC flag in the EVTS register. The query form *OPC? returns 1 when complete, but does not affect the EVTS register.
Example	*OPC? <i>1</i>
Related commands	

PRELIMINARY INFORMATION

CONS (Console Mode)

Group	Interface commands
Action	Set/Query
Syntax	CONS(?) { <i>m</i> }
Description	Set (query) the Console mode {to <i>m</i> }. CONS 1 causes each character received to be returned to the host computer.
Allowed range	$m \in \{0 \text{ (disabled)}, 1 \text{ (enabled)}\}$
Reset (*RST) value	0
Power-on value	0
Example	CONS 1 <i>1</i>
Related commands	

PRELIMINARY INFORMATION

***IDN (Identify)**

Group	Interface commands
Action	Query only
Syntax	*IDN?
Description	<p>Read the device identification string. This string is formatted as:</p> <p style="padding-left: 40px;">Signals and Systems for Physics, model SK305, hw Rppx, fw Rqqy, s/n dddddd.</p> <p>In this string, SK305 is the model number, Rnnx and Rppy are revision numbers identifying the hardware or the firmware versions and dddddd is the 6-digit serial number.</p>
Example	<p>*IDN?</p> <p><i>Signals and Systems for Physics, model SK305, hw R24B, fw R24A, s/n 123456.</i></p>
Related commands	

PRELIMINARY INFORMATION

LINS (Last Instrument Error Status)

Group	Status reporting commands
Action	Query only
Syntax	LINS?
Description	Query the last execution instrument error. LINS? returns the unique code number associated with this event.
Valid codes are	<p>0 \longleftrightarrow no execution error since last LINS?;</p> <p>1 \longleftrightarrow on-chip ADC error;</p> <p>10 \longleftrightarrow detected hardware is in invalid condition;</p> <p>20 \longleftrightarrow some parameters have been be adapted or clamped;</p> <p>21 \longleftrightarrow some functionalities have been be disabled.</p>
Power-on value	0
Example	LINS? 0
Related commands	LCMD, LEXE, LURQ.

PRELIMINARY INFORMATION

LURQ (Last User Request Status)

Group	Interface commands
Action	Query only
Syntax	LURQ?
Description	Query the last user request. LURQ? returns the unique code number associated with this event.
Valid codes are	0 \leftrightarrow No User request since last LURQ?
Power-on value	0
Example	LURQ? 0
Related commands	LCMD, LEXE, LINS.

PRELIMINARY INFORMATION

LCMD (Last Command Error Status)

Group	Interface commands
Action	Query only
Syntax	LCMD?
Description	Query the last command error. LCMD? returns the unique code number associated with this error.
Valid codes are	<p>0 ↔ no execution error since last LCMD?</p> <p>1 ↔ illegal (unknown) command.</p> <p>2 ↔ illegal query.</p> <p>3 ↔ illegal set (read-only command).</p> <p>4 ↔ extra parameter.</p> <p>5 ↔ missing parameter.</p> <p>6 ↔ null command.</p>
Power-on value	0
Example	*RST?;LCMD? 2
Related commands	LURQ, LEXE, LINS.

PRELIMINARY INFORMATION

LEXE (Last Execution Error Status)

Group	Interface commands
Action	Query only
Syntax	LEXE?
Description	Query the last execution error. LEXE? returns the unique code number associated with this error.
Valid codes are	<p>0 ↔ no execution error since last LEXE?</p> <p>1 ↔ invalid parameter.</p> <p>2 ↔ argument value out-of-range.</p> <p>3 ↔ some parameters have been adapted or clamped.</p> <p>4 ↔ a conflict due to the current operation has been avoided.</p> <p>5 ↔ no change upon executing the command.</p> <p>6 ↔ the operation was aborted due to a fault condition.</p>
Power-on value	0
Example	<pre>CONS2;LEXE?;LEXE?</pre> <p>1</p> <p>0</p>
Related commands	LURQ, LCMD, LINS.

TERM (Response Termination)

Group	Interface commands
Action	Set/Query
Syntax	TERM(?) { <i>m</i> }
Description	<p>Set (query) the termination sequence {to <i>m</i>}.</p> <p>The termination sequence is appended to all query responses sent by the instrument. It is constructed of ASCII character(s) <CR> (carriage return) or <LF> (line feed).</p>
Allowed range	<p>$m \in \{1, 2, 3, 4\}$ where :</p> <p>1 \longleftrightarrow <CR> character appended, 2 \longleftrightarrow <LF> character appended, 3 \longleftrightarrow both <CR> and <LF> characters appended, 4 \longleftrightarrow no character appended.</p>
Power-on value	3
Reset (*RST) value	3
Example	TERM? 3
Related commands	

2.6 Memory commands

The Memory commands allow the User to save and recall the instrument's settings in non-volatile memory.

*RCL (Recall Settings)

Group	Memory commands
Action	Set only
Syntax	*RCL
Description	Recall the settings stored in the non-volatile memory.
Example	*RCL
Related commands	*RST, *SAV.

PRELIMINARY INFORMATION

***SAV (Save Current Settings)**

Group	Memory commands
Action	Set only
Syntax	*SAV
Description	Save the current settings in the non-volatile memory.
Example	*SAV
Related commands	*RCL, *RST.

PRELIMINARY INFORMATION

3 Status Model

The complete block diagram of the status register array is available online at the related product page. There are four categories of registers in this model :

Last Event registers These four read registers (LINS, LCMD, LURQ and LEXE) store the last event that they monitor. A query command i) return the last registered event since the previous query and ii) clears the register's content.

Condition registers These read-only registers correspond to the real-time condition of some underlying physical properties under monitoring. Queries return the latest value of the property, and have no other effect.

Condition register names end with C.

Status registers These read-only registers record the occurrence of defined events. If the event occurs, the corresponding status bit is set to 1. Upon querying a status register, any set bits within it are cleared. These are sometimes known as sticky bits since once set, a bit can only be cleared by reading its value. Status register names end with S.

Enable registers These read/write registers define a bitwise mask for their corresponding status register. If any bit position is set in a status register while the same bit position is also set in the enable register, then the corresponding summary bit is set in either the Event Summary or Master Summary register. Enable register names end with E.

PRELIMINARY INFORMATION

3.1 Master Summary Status (MSTS)

The Master Summary Status (MSTS) is the top-level summary register of the status model. When masked by the Master Summary Status Enable (MSTE) register, a bit set in the Status Byte causes the /STATUS signal to be asserted on the DIN41612 connector. This register is queried with the `MSTS?[n]` command.

Weight $n = 2^i$	Bit i	Flag	Description
128	7	OVL	Overload Summary Bit. Indicates whether one or more of the enabled flags in the Overload Status register is true.
64	6	INS	Instrument Summary Bit. Indicates whether one or more of the enabled flags in the Instrument Status register is true.
32	5	RFU	Undefined (read 0).
16	4	RFU	Undefined (read 0).
8	3	RFU	Undefined (read 0).
4	2	EVT	Event Summary Bit. Indicates whether one or more of the enabled flags in the Event Status register is true.
2	1	COM	Communication Summary Bit. Indicates whether one or more of the enabled flags in the Communication Status register has become true.
1	0	MSS	Master Summary Status. Indicates whether one or more of the enabled status messages in the Status Byte register is true.

3.2 Master Summary Enable (MSTE)

Each bit in the MSTE register corresponds one-to-one with a bit in the MSTS register, and acts as a bitwise AND of the MSTS flags to generate the MSS flag. Bit 0 of the MSTE is undefined—setting it has no effect, and reading it always returns 0. This register is set and queried with the `MSTE(?)` command and cleared at power-on.

3.3 Event Status (EVTS)

The Event Status register consists of 8 event flags. These flags are set by the corresponding event, and cleared only by reading or with the *CLS command ("sticky bits"). Querying the single bit i with the command EVTS? n where the bit-mask $n = 2^i$ will only clear the bit i . For instance, issuing the command EVTS?128 will clear the bit 7 (INS) only.

Weight $n = 2^i$	Bit i	Flag	Description
128	7	INS	Instrument event. Indicates whether one or more of the enabled flags in the Instrument Status register is true.
64	6	URQ	User Request event. Indicates that a User request has occurred. The request code can be queried with LURQ?
32	5	TXQ	Transmission Buffer event. Indicates that the TX buffer has been cleared.
16	4	RXQ	Reception Buffer event. Indicates that the RX buffer has been cleared.
8	3	EXE	Execution Error event. Indicates an error in a command that was successfully parsed. The error code can be queried with LEXE?
4	2	CMD	Command Error event. Indicates an error detected by the command parser. The error code can be queried with LCMD?
2	1	OPC	Operation Complete. Set by the *OPC command.
1	0	PON	Power On event. Indicates that an off-to-on transition has occurred.

3.4 Event Enable (EVTE)

Each bit in the EVTE register corresponds one-to-one with a bit in the EVTS register, and acts as a bitwise AND of the EVTS flags to generate the EVT flag in the Master Summary Status (MSTS) register. This register is set and queried with the EVTE command and cleared at power-on. For instance, issuing the command EVTE 128 enable the bit 7 (INS) only.

3.5 Instrument Status (INSS)

The Instrument Status register consists of 8 event flags. These flags are set by the corresponding event, and cleared only by reading or with the *CLS command ("sticky bits"). Querying the single bit i with the command INSS? n where the bit-mask $n = 2^i$ will only clear the bit i . For instance, issuing the command INSS?1 will clear the bit 0 only.

Weight $n = 2^i$	Bit i	Flag	Description
128	7	RFU	Undefined (read 0).
64	6	RFU	Undefined (read 0).
32	5	RFU	Undefined (read 0).
16	4	TPO	The current source has tripped off due to a fault detection (bit set).
8	3	OPN	Open-circuit detected at output (bit set).
4	2	ENA	The TEC current source is on (bit set).
2	1	IKS	Internal 10-MHz clock source used. The module is not synchronized to the platform's timebase. <i>The platform's timebase synchronization feature is not yet implemented. The bit is therefore always set (read 1).</i>
1	0	PUV	At least, one power supply is under its low-level threshold (bit set).

3.6 Instrument Enable (INSE)

Each bit in the INSE register corresponds one-to-one with a bit in the INSS register, and acts as a bitwise AND of the INSS flags to generate the INS flag in the Master Summary Status (MSTS) register. This register is set and queried with the INSE command and cleared at power-on.

3.7 Instrument Condition (INSC)

Each bit in the INSC register corresponds one-to-one with a bit in the INSS register. The bits in the INSC register reflect the real-time values of their corresponding signals. Reading the entire register, or individual bits within it, does not affect the value of INSC. This register is queried with the INSC command and cleared at power-on.

3.8 Overload Status (OVLS)

The Overload Status register consists of 8 event flags. These event flags are set by the corresponding event, and cleared only by reading or with the *CLS command ("sticky bits"). Querying the single bit i with the command OVLS? n where the bit-mask $n = 2^i$ will only clear the bit i . For instance, issuing the command OVLS?2 will clear the bit 1 only.

Weight $n = 2^i$	Bit i	Flag	Description
128	7	RFU	Undefined (read 0).
64	6	RFU	Undefined (read 0).
32	5	RFU	Undefined (read 0).
16	4	OVT	The die temperature of the power stage is excessive (bit set).
8	3	VTN	The output voltage reaches its lower limit (bit set).
4	2	VTP	The output voltage reaches its upper limit (bit set).
2	1	ILN	The operating current is below the negative limit (bit set).
1	0	ILP	The operating current is above the positive limit (bit set).

3.9 Overload Enable (OVLE)

Each bit in the OVLE register corresponds one-to-one with a bit in the OVLS register, and acts as a bitwise AND of the OVLS flags to generate the OVL flag in the Master Summary Status (MSTS) register.

3.10 Overload Condition (OVLC)

Each bit in the OVLC register corresponds one-to-one with a bit in the OVLS register. The bits in the OVLC register reflect the real-time values of their corresponding signals. Reading the entire register, or individual bits within it, does not affect the value of OVLC. This register is queried with the OVLC command and cleared at power-on.

3.11 Communication Status (COMS)

The Communication Status register consists of 8 event flags. These flags are set by the corresponding event, and cleared only by reading or with the *CLS command ("sticky bits"). Querying the single bit i with the command COMS? n where the bit-mask $n = 2^i$ will only clear the bit i .

Because the COMS register is not used in the SK305, querying this register always returns 0. Therefore, the corresponding summary bit in the MSTS register (bit COM) is never set whatever the value of the COME register.

Weight $n = 2^i$	Bit i	Flag	Description
128	7	RFU	Undefined (read 0).
64	6	RFU	Undefined (read 0).
32	5	RFU	Undefined (read 0).
16	4	RFU	Undefined (read 0).
8	3	RFU	Undefined (read 0).
4	2	RFU	Undefined (read 0).
2	1	COL	Bus collision.
1	0	PRY	Parity violation.

3.12 Communication Enable (COME)

Each bit in the COME register corresponds one-to-one with a bit in the COMS register, and acts as a bitwise AND of the COMS flags to generate the COM flag in the Master Summary Status (MSTS) register. This register is set and queried with the COME command and cleared at power-on.

PRELIMINARY INFORMATION

3.13 Last Command Error (LCMD)

The LCMD register holds the last error detected by the command parser. The related error code can be retrieved by the command LCMD?. When such an error is detected, the corresponding bit in the Event Status register is set (bit CMD in EVTS).

3.14 Last Execution Error (LEXE)

The LEXE register holds the last error detected during the execution of a command. The related error code can be retrieved by the command LEXE?. When such an error is detected, the corresponding bit in the Event Status register is set (bit EXE in EVTS).

3.15 Last Instrument Error (LINS)

The LINS register holds the last error detected during the operation of the instrument. The related error code can be retrieved by the command LINS?. When such an error is detected, the corresponding bit in the Event Status register is set (bit INS in EVTS).

3.16 Last User Request (LURQ)

The LURQ register holds the last User's request. The related request code can be retrieved by the command LURQ?. When such a request is reported, the corresponding bit in the Event Status register is set (bit URQ in EVTS).

Because the LURQ register is not used in the SK305, querying this register always returns 0 and the corresponding summary bit in the Event Status register is never set (bit URQ in EVTS).

4 Index of commands

Instrument Configuration commands

- EXTE (External Control State), 15
- FFWE (External Current Control State), 16
- ITPO (Trip-off upon Current Limiting), 18
- MANE (Manual Control State), 14
- MONS (MONO Output Selection), 20
- TECE (TEC Output State), 17
- VTPO (Trip-off upon Overvoltage Detection), 19

Instrument Monitoring commands

- RMON (Remote Monitoring), 21
- STME (Data Streaming Enable), 24
- STMN (Number of Streamed Measurements), 25
- STMS (Streamed Channels Selection), 23
- TDIE (Die Temperature), 22

Instrument Settings commands

- FFWG (Feedforward Gain), 13
- ILMN (Negative Current Limit), 10
- ILMP (Positive Current Limit), 9
- MANS (Manual Current Control), 8
- VTHN (Negative Voltage Threshold), 12
- VTHP (Positive Voltage Threshold), 11

Interface commands

- *IDN (Identify), 42

- *OPC (Operation Complete), 40

- *RST (Reset), 39

- CONS (Console Mode), 41

- LCMD (Last Command Error Status), 45

- LEXE (Last Execution Error Status), 46

- LINS (Last Instrument Error Status), 43

- LURQ (Last User Request Status), 44

- TERM (Response Termination), 47

Memory commands

- *RCL (Recall Settings), 48

- *SAV (Save Current Settings), 49

Status Reporting commands

- *CLS (Clear Status Registers), 26

- COME (Communications Enable), 32

- COMS (Communications Status), 31

- EVTE (Event Enable), 30

- EVT S (Event Status), 29

- INSC (Instrument Condition), 38

- INSE (Instrument Enable), 37

- INSS (Instrument Status), 36

- MSTE (Master Summary Enable), 28

- MSTS (Master Summary Status), 27

- OVLC (Overload Condition), 35

- OVLE (Overload Enable), 34

- OVLS (Overload Status), 33

5 Document Revision History

5.1 Version Number

This document is identified by SK305-SU01-P25A.

5.2 Revision History

P25A (2025-01-07)

Initial version.

PRELIMINARY INFORMATION